HUMANTECHNIK

iBell

D

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Seite 2

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GB

Operation Instructions

Page 12

Call detector for mobile phones

F

Mode d'emploi

Page 22

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Congratulations on purchasing your »iBell«. We hope it will provide you with many years of service. You have chosen a modern and reliable system. Please read through these operating instruction carefully to be able to start the unit correctly and to familiarise yourself with all of the systems features.

Standard components

Please check if all following components are included:

- »iBell«
- 3x 1,5 V AA batteries
- Connection cable
- Operating Instructions
- Warranty card

If any parts are missing please immediately contact your dealer or the manufacturer directly.

Safety instructions:


- Place the device on a non-slip surface that does not vibrate.
- Ensure the leads cannot cause any accidents.
- Keep the device away from the following:
 - Sources of heat
 - Direct sunlight
 - Humidity
- Protect your device from moisture, dust, aggressive liquids and vapours.
- Only connect the power leads to the appropriate jacks.
- Only connect authorised accessories.
- Never:
 - open the device yourself;
 - carry the device by its leads.

Table of contents

Standard components	Page 12
Safety instructions	Page 12
Controls	Page 14, 15
How it works	Page 16
Starting device	Page 16
Inserting the batteries	Page 16
Ready indicator/ battery indicator	Page 16
Positioning the »iBell«	Page 17
Positioning the mobile phone	Page 17
Call/text message received confirmation	Page 17
Connecting the »iBell« with the RF <i>lisa</i> signalling system	Page 18
Sensitivity settings for the acoustic RF telephone transmitter	Page 18
Connecting additional modules to the »iBell«	Page 18
Troubleshooting	Page 19
Maintenance and care	Page 20
Warranty	Page 20
Environmental note/ disposal	Page 21
Technical data	Page 21
CE symbol	Page 21

Controls

Ready indicator/ battery indicator

LED panel on reverse

Battery compartment

Controls

Component	Function
Battery compartment	Contains 3x AA batteries.
Ready indicator/ battery indicator	Shows that the »iBell« is on.
Reset button	Confirms a call or text message has been received.
LED panel on front	Blinks as soon as the mobile phone has received a call or a text message.
LED panel on reverse	Blinks as soon as the mobile phone has received a call or a text message.
Connection jack	To connect the »iBell« with the <i>lisa</i> signalling system or additional modules such as a vibrating pillow.

How it works

When the mobile phone in the »iBell« receives a text message or call, the »iBell« uses vibration sensors to detect the mobile phone is vibrating and triggers a blue chaser light to attract attention.

For the first 2 minutes after a call is received, the light blinks continually. After 2 minutes, blinking slowly decreases and eventually disappears completely.

Starting device

Inserting the batteries

The device needs three 1.5 V batteries (AA type) to work.

To insert the batteries, proceed as follows:

1. Press down the lid of the battery compartment on the back of the »iBell«.
2. Now insert three AA type batteries into the battery compartment.
Ensure the right polarity is selected for the batteries.
3. Now close the battery compartment lid.

As soon as the batteries have been inserted, a signal will sound and a blue chaser light on the »iBell« will briefly blink.

Ready indicator/battery indicator:

Ready indicator/battery indicator:

The blue LED at the top left will indicate that the device is on.

- If it blinks every 15 seconds it means that:
your batteries are full.
- Constant blinking means that:
the batteries are weak, please exchange them.

If the batteries are too weak, the device will switch off automatically.

Positioning the »iBell«

Place the »iBell« on a flat, non-vibrating surface.

Suitable surfaces: tables, shelves etc.

Unsuitable surfaces: tumble dryers, washing machines etc.

Positioning the mobile phone

Positioning the mobile phone

1. Switch your mobile phone to vibration alarm.
If your mobile phone has different vibration rhythms, select the alarm with the longest and strongest vibration.
2. Carefully place your mobile phone into the »iBell's« housing.
Please note: if the mobile phone is not inserted carefully, the vibration alarm could make the »iBell« start to blink.
If that happens, press the reset button to acknowledge the signal.

Call/text message received confirmation


After receiving a call or a text message, you must confirm it has been received so that the »iBell« can report new calls or text messages.

Press the reset button on the front of the »iBell« to acknowledge the signal. The device has now been reset and can report new calls and text messages again.

Connecting the »iBell« with the RF *lisa* signalling system

The »iBell« can be connected to the RF *lisa* signalling system.

A signalling system consists of at least one transmitter and one receiver. The acoustic RF telephone transmitter (not included) can be connected to the »iBell« with the connecting lead supplied.


The acoustic RF telephone transmitter can now accept the »iBell's« signals and change them into radio pulses. These radio pulses (868.35 MHz) are sent to the receiver by wireless transmission.

Sensitivity settings for the acoustic RF telephone transmitter

Set the sensitivity control on the bottom of the acoustic RF telephone transmitter to the highest level (turn the switch to the right until it goes no further). The transmitter will now react to the »iBell's« signals.

Connecting additional modules to the »iBell«

Signal emitter modules (vibrating pillow, flash module MF-1, acoustic module MA-1 or switch module MS-1) can be connected with the »iBell« via its power supply jack.

Plug one end of the supplied lead in the »iBell's« connection jack and the other end of the lead in the module's connection jack.

Troubleshooting

Symptoms	Possible causes/remedy
The »iBell« does not blink when a call/text message has been received.	Check whether the batteries have been inserted correctly and are fully charged. Check whether the longest and strongest vibration alarm has been set on your mobile phone.
»iBell« blinks although no call/text message has been received.	The device is very close to a strong wireless signal. Press the reset button and put the mobile phone back in the housing. The device was placed on an uneven or vibrating surface (e.g. worktop with kitchen equipment). The mobile phone was not placed in the housing carefully.
The blue LED at the top left is constantly blinking.	Battery indicator: insert new batteries.
I have connected the »iBell« with the wireless <i>lisa</i> signalling system, but it is not triggering.	Ensure that all the leads are connected properly. Contact your dealer.

Maintenance and care

If the unit does become dirty, simply wipe it clean with a soft, damp cloth. Never use spirits, thinners or other organic solvents. Do not set up the »iBell« where it will be exposed to full sunlight for long periods. In addition, it must be protected against excessive heat, moisture and severe mechanical shocks.

Note: This product is **not** protected against splash water. Do not place any containers filled with water, such as flower vases, or anything with an open flame, such as a lit candle, on or near the product. Please make sure that the batteries are not exposed to excessive heat from such as sunlight, fire or anything similar.

Warranty

This product is very safe to operate and provides highly reliable service. Should a malfunction occur despite the unit having been set up and operated correctly, please contact your dealer or the manufacturer directly.

This warranty covers the repair of the product and returning it to you free of charge. It is essential that you send in the product in its original packaging, so do not throw the packaging away. The warranty does not apply to damage caused by incorrect handling or attempts to repair the unit by people not authorised to do so (destruction of the seal on the unit). Repairs will only be carried out under warranty if the completed warranty card is returned accompanied by a copy of the dealer's invoice/till receipt.

Always specify the product number in any event.

Environmental note / disposal

Disposal of used electrical and electronic equipment (applicable in the countries of the European Union and in other European countries which operate a separate collection system for these devices).

The symbol on the product or its packaging indicates that this product should not be treated as normal household waste. It should be handed in to a collection center for recycling electrical and electronic equipment. By playing your part in ensuring the correct disposal of this product, you will be safeguarding the environment and the health of other people. The environment and health are jeopardized by incorrect disposal. Material recycling helps to reduce the consumption of raw materials. You can obtain further information about recycling this product from your local council, from local waste disposal operators and from the shop where you bought the product.

Technical data

Power supply:	3x 1,5 V AA batteries
Height:	102 mm
Width:	93 mm
Depth:	103 mm
Weight:	224 g incl. batteries

CE symbol

This device satisfies the following EU directives:

- 2002/95/EG RoHS directive
- 2002/96/EG WEEE directive
- 1999/5/EG R&TTE directive
- 2004/108/EG EMC directive
- 2006/95/EG Low voltage directive

Compliance with the directives listed above is confirmed by the CE seal on the device. CE compliance declarations are available on the Internet at www.humantechnik.com.

Technical specifications subject to change without notice.

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Humantechnik GHL AG